

ET 350

Changes of state in the refrigeration circuit



Description

- refrigeration circuit demonstrated clearly
- transparent components offer insights into the changes of state
- energetic analyses of the refrigeration cycle
- value of Global Warming Potential $GWP=1$

In a compression refrigeration system a refrigerant flows through the refrigeration circuit and is subject to different changes of state. Here, the physical effect is used that during the transition of the refrigerant from a liquid to a gaseous state energy is required which is removed from the environment (evaporation enthalpy).

The experimental unit ET 350 represents a typical refrigeration circuit consisting of a hermetic piston compressor, condenser, expansion valve and evaporator. The evaporator and condenser are transparent to provide good monitoring of the phase transition process during evaporation and condensation. The operation of the float valve as expansion valve is also easy to observe. Before the entry into the evaporator the aggregate state of the refrigerant can be monitored at a sight glass. A water circuit cools the condenser or supplies the cooling load for the evaporator.

Cold and hot water and refrigerant flows are adjustable. The low pressure level of the refrigerant used permits the use of an evaporator and condenser out of glass.


Temperatures and pressures are recorded and displayed. The key points of the cyclic process can be read and entered into a log p-h diagram. The power of the compressor and flow rates of the water flows and the refrigerant are also indicated.

Learning objectives/experiments


- design and operation of a compression refrigeration system
- observe the evaporation and condensation of the refrigerant
- represent and understand the refrigeration cycle in the log p-h diagram
- energy balances
- calculation of the coefficient of performance

ET 350


Changes of state in the refrigeration circuit



1 pressure switch, 2 flow meter, 3 condenser, 4 expansion valve, 5 compressor, 6 evaporator, 7 sight glass, 8 temperature display, 9 manometer



1 condenser, 2 expansion valve, 3 compressor, 4 evaporator; T temperature, P pressure, E electrical power, F flow rate, PSH pressure switch; blue: low pressure, red: high pressure, green: water



Expansion valve in the shape of a float valve: 1 float lifts the needle from the valve seat, 2 needle, 3 refrigerant escapes, 4 tank

Specification

- [1] demonstration of the processes in a refrigeration circuit
- [2] for better process monitoring the evaporator and condenser are of transparent design
- [3] evaporator and condenser with pipe coil
- [4] expansion valve in the shape of a float valve
- [5] pressure switch to protect the compressor
- [6] temperature sensor, power meter, manometer in the refrigeration circuit, flow meter for hot and cold water and refrigerant
- [7] water supply via the lab network or via WL 110.20 Water chiller to ensure a maximum water temperature of 16°C
- [8] safety valves at the evaporator and condenser
- [9] refrigerant R1233zd, GWP: 1

Technical data

Hermetic piston compressor
 ■ capacity: 18,3cm³

Evaporator capacity: approx. 2800mL
 Condenser capacity: approx. 2800mL

Refrigerant
 ■ R1233zd
 ■ GWP: 1
 ■ filling volume: 1,2kg
 ■ CO₂-equivalent: 0t

Measuring ranges
 ■ temperature: 8x -20...200°C
 ■ pressure: 2x -1...1,5bar
 ■ flow rate:
 ▶ 2x 0...48L/h (water)
 ▶ 0...700L/h (refrigerant)
 ■ power: 0...1200W

230V, 50Hz, 1 phase
 230V, 60Hz, 1 phase
 120V, 60Hz, 1 phase
 UL/CSA optional
 LxWxH: 1200x500x900mm
 Weight: approx. 110kg

Required for operation

water connection (min. 48L/h, water temperature max. 16°C), drain or WL 110.20

Scope of delivery

- 1 experimental unit
- 1 set of hoses
- 1 set of instructional material

ET 350

Changes of state in the refrigeration circuit

Optional accessories

| | | |
|-----------|-----------|--------------------|
| 020.30009 | WP 300.09 | Laboratory trolley |
| 060.11020 | WL 110.20 | Water chiller |